

Unit Dose Presoak

Enzymatic formula breaks down protein and starch films



Effective and efficient

- Presoaking flatware results in consistently clean and spotless results with every wash cycle
- Saves time and labor by helping to reduce need for rewash

Safe to use

- Safer product handling
- Requires no PPE

Easy on the environment

- Concentrated solid tablets contribute to reduced water, energy, and plastic waste

Unit Dose tablets weigh less than standard ready-to-use solutions



DIRECTIONS FOR USE:

1. Tear open a packet and drop the tablet into a presoak tub. Wash hands after handling product.
2. Fill up the tub with warm water (110° – 120° F) and allow the tablet to dissolve.
3. Place dirty flatware into the solution and allow to soak.

Unit Dose Presoak	
PACK SIZE	48 ct
DILUTION RATIO	1 tablet/ tub
PICK CODE	6102700
SUPC CODE	7303533



We thought small so you can think big.
 Ready to learn more about our full tablet cleaning program?
 Visit cleanwithkeystone.com